

Entry points: Python API · CLI · Workbench · strict BL19B2 batch workflow

Inputs

2D detector frames
and instrument headers

External 1D profiles
and correction metadata

NIST SRM 3600,
water, or supplied
reference curve

Shared processing

Parse, normalize, and
validate inputs

Reduce detector data and
integrate with pyFAI

Estimate K and propagate
supported uncertainty

Check intensity state
and correction history
before calibrated output

Outputs

Absolute 1D profiles
and calibrated 2D data

canSAS XML and
NXcanSAS HDF5

Calibration records,
source hashes, QC,
and run reports